



Test Report

No. CANEC0901624701

Date: 16 Apr 2009

Page 1 of 3

Distributor Name: LORIKA CZ s.r.o., Holická 312, 77200 Olomouc, Czech Republic

The following sample(s) was/were submitted and identified on behalf of the clients as :

Polyester swab

SGS Job No. : 10935466 - SZ

SGS Internal Reference No. : 2.1

Date of Sample Received : 11 Apr 2009

Testing Period : 11 Apr 2009 - 16 Apr 2009

Test Requested : Selected test(s) as requested by client.

Test Method : Please refer to next page(s).

Test Results : Please refer to next page(s).

Conclusion : Based on the performed tests on submitted sample(s), the results **comply with the RoHS Directive 2002/95/EC and its subsequent amendments.**Signed for and on behalf of
SGS-CSTC Ltd.

Huang Fang, Sunny
Sr. Engineer

This document is issued by the Company under its General Conditions of Service printed overleaf or available on request and accessible at http://www.sgs.com/terms_and_conditions. Attention is drawn to the limitation of liability, indemnification and jurisdiction issues defined therein. Any holder of this document is advised that information contained herein reflects the Company's findings at the time of its intervention only and within the limits of Client's instructions, if any. The Company's sole responsibility is to its Client and this document does not generate parties to a transaction from exercising all their rights and obligations under the transaction documents. Any unauthorized alteration, forgery or falsification of the content or appearance of this document is unlawful and offenders may be prosecuted to the fullest extent of the law. Unless otherwise stated, the results shown in this test report refer only to the sample(s) tested and such sample(s) are retained for 30 days only. This document cannot be reproduced except in full, without prior approval of the Company.



SGS-CSTC (Guangdong) Technical Service Co., Ltd. 198 Kexu Road, SCENTECH Park, Guangzhou Economic & Technology Development District, Guangzhou, China 510663 t: (86-20) 82155555 f: (86-20) 82075125 www.cn.sgs.com
Guangzhou SGS-CSTC Technical Service Co., Ltd. 中国·广州·经济技术开发区科学城科珠路198号 邮编: 510663 t: (86-20) 82155555 f: (86-20) 82075125 e: sgs.china@sgs.com

Member of the SGS Group (SGS SA)



ANALABS TECHNOLOGIES (WUXI) CO.LTD.

No.789 Zhongnan Road, Block A, 5th Floor Yangming
Digital Electronic Park Wuxi, Jiangsu Province P.R.China, P.C.214024
(Tel):86 510-85430236 (Fax):86 510-85430235

REPORT

Distributor Name: LORIKA CZ s.r.o., Holická 312, 77200 Olomouc, Czech Republic

Sample Description : One (1) sample of Foam Swabs

The sample consisted of five packages of Foam Swabs and marked :

ini ompressed

Foam Swab

P/N: CM-FS750

Kleentips Swab 750

PKG: 100pcs/ Bag 500pcs/pack

On analysis in accordance with Seagate's Doc.No.35344 001 Rev.P Method, the following results were obtained:

Test	Result	Seagate Doc.No. 35348-001 Rev.C Specification
Chloride (Cl), ug/tip	0.15	≤ 0.90
Sulphate (So ₄), ug/tip	0.02	≤ 0.40
Total Anions, ug/tip	0.17	≤ 1.3
Non-Volatile Residue, NVR, mg/swab (20ml D.I. Water Extraction, at ambient temperature for 24 Hours)	< 0.01	≤ 0.15
Non-Volatile Residue, NVR, mg/swab (20ml IPA, at ambient temperature for 24 Hours)	0.03	≤ 0.40

FT-IR/ATR Analysis of the
Non-Volatile Residue Extract,
(Hexane)Amide, DOP and
Silicone Oil were not detected
(Refer to Spectrum attached)No detectable of Phthalates
and Amides.
(< 1.3ng/tip of silicone oil)LPC Test (≥ 0.5um)
Particle Counts/swab tip

Data #1	Data #2	Data #3	Average	
3,911	4,789	3,112	3,937	< 25,000

Remarks : (1). "<" - Less than
(2). "≤" - Smaller than or equal to
(3). "≥" - Greater than or equal to

SIM KOK HAI
DIRECTOR
ANALABS TECHNOLOGIES (WUXI) CO.LTD.

LORIKA CZ s.r.o.

Holická 312

779 00 Olomouc

Tel. +420 581 112 410

Fax +420 581 112 466

e-mail: lorika@lorika.czwww.lorika.cz

LORIKA CZ s.r.o., Náchod

Tel./fax +420 491 423 618

e-mail: nachod@lorika.cz

LORIKA CZ s.r.o., Teplice

Tel./fax +420 417 577 497

e-mail: teplce@lorika.cz

LORIKA SLOVAKIA s.r.o.

Langsfeldova 2

036 01 Martin

Tel. +421 434 220 057

Fax +421 434 287 277

e-mail: lorika@lorika.skwww.lorika.sk

ARISO® SWABS

Physical Characteristics:

To avoid the static electricity to damage and high Class Cleanroom demand, Please choose Cleanmo Anti-static and Cleanroom swabs.

P/N:	Head Material	Head width	Head Thickness	Head length	Handle Material	Handle width	Handle Length	Total length	Packaging
CM-FS750 (TX750B/D)	Sponge	3.2mm	3.2mm	10mm	Polypropylene	2.5mm	66mm	76mm	500swabs/bag 5 inner bags of 100 swabs 50bags/ctn
CM-FS757 (TX757B/E)	Sponge	3.2mm	2.6mm	10mm	Polypropylene	2.2mm	60mm	70mm	500swabs/bag 5 inner bags of 100 swabs 50bags/ctn
CM-FS742 (TX742B/C)	Sponge	3.6mm	3.5mm	12mm	Polypropylene	3.0mm	59mm	71mm	500swabs/bag 5 inner bags of 100 swabs 50bags/ctn
CM-FS741 (TX741B)	Sponge	3.2mm	3.0mm	11mm	Polypropylene	2.2mm	57mm	68mm	500swabs/bag 5 inner bags of 100 swabs 50bags/ctn
CM-FS751 (TX751B/S)	Sponge	4.0mm	4.0mm	15mm	Polypropylene	2.4mm	51mm	66mm	500swabs/bag 5 inner bags of 100 swabs 50bags/ctn
CM-FS740 (TX740B)	Sponge	6.8mm	5.0mm	18mm	Polypropylene	3.0mm	145mm	163mm	500swabs/bag 5 inner bags of 100 swabs 10bags/ctn
CM-FS746 (TX746S)	Sponge	6.8mm	5.0mm	18mm	Polypropylene	3.0mm	74mm	92mm	500swabs/bag 5 inner bags of 100 swabs 20bags/ctn
CM-FS712 (TX712A)	Sponge	11mm	7.0mm	25mm	Polypropylene	5.2mm	103mm	128mm	50swabs/bag 100bags/ctn or 50bags/ctn
CM-FS706 (TX706A)	Sponge	11mm	7.0mm	27mm	Polypropylene	6.3mm	54mm	81mm	50swabs/bag 50bags/ctn
CM-FS707 (TX707A)	Sponge	16mm	8.0mm	25mm	Polypropylene	6.6mm	106mm	131mm	50swabs/bag 50bags/ctn
CM-FS718	Sponge	13mm	6.0mm	25mm	Polypropylene	5.0mm	75mm	128mm	50swabs/bag 50bags/ctn

P/N:	Head Material	Head width	Head Thickness	Head length	Handle Material	Handle width	Handle Length	Total length	Packaging
CM-FS708 (TX708A)	Sponge	21mm	8.0mm	26mm	Polypropylene	6.6mm	106mm	132mm	50swabs/bag 50bags/ctn
CM-FS920	Sponge	6.0mm	5.0mm	11mm	Polypropylene	4.0mm	57mm	78mm	100swabs/bag 100bags/ctn
CM-FS719	Sponge	3.0mm	3.0mm	11mm	Polypropylene	2.5mm	57mm	68mm	500swabs/bag 5*100swabs/bag 50bags/ctn
CM-PS743 (TX743B)	Double-Layer Polyester	4.0mm	3.6mm	12mm	Polypropylene	3.0mm	58mm	70mm	500swabs/bag 5*100swabs/bag 50bags/ctn
CM-PS758L (TX758L)	Single-Layer Polyester	3.2mm	2.2mm	10mm	Polypropylene	2.2mm	60mm	70mm	500swabs/bag 5*100swabs/bag 50bags/ctn
CM-PS758LM	Single-Layer Microfiber	3.2mm	2.2mm	10mm	Polypropylene	2.2mm	60mm	70mm	500swabs/bag 5*100swabs/bag 50bags/ctn
CM-PS759 (TX759B)	Microfiber	2.2mm	1.2mm	10mm	Polypropylene	2.2mm	62mm	72mm	500swabs/bag 5*100swabs/bag 50bags/ctn
CM-PS714 (TX714A)	Double-Layer Polyester	13mm	3.0mm	24mm	Polypropylene	5.2mm	104mm	128mm	100swabs/bag 50bags/ctn
CM-PS761 (Tx761)	Double-Layer Polyester	7.5mm	3.0mm	16mm	Polypropylene	3.2mm	148mm	164mm	500swabs/bag 5*100swabs/bag 10bags/ctn
CM-PS766	Double-Layer Polyester	7.5mm	3.0mm	16mm	Polypropylene	3.2mm	72mm	92mm	500swabs/bag 5*100swabs/bag 20bags/ctn
CM-PS765 (Tx765)	Singel-Layer Polyester	3.0mm	2.5mm	9mm	Polypropylene	3.2mm	143mm	153mm	500swabs/bag 5*100swabs/bag 10bags/ctn
CM-PS707M	Microfiber	13.0mm	3.0mm	25mm	Polypropylene	6.6mm	105mm	130mm	100swabs/bag 50bags/ctn